

LOLER Report

<p style="text-align: right;">C1027</p> <p>AIG Control Systems Ltd Otter Pool Barn West Stowford Barton Barn, West Stowford, Ilfracombe Devon, England EX34 8PJ</p> <p>07970 538998</p>	<p style="text-align: center;">Visit Details</p> <p>Visit Type: LOLER Order No: 94541/ Date/Time: 14/11/2025 11:29 Time Completed: 11:53</p>
<p>Job Details</p> <p>LOLER inspection visit due, 2025</p>	<p style="text-align: center;">Lift Details</p> <p>Asset /Location: Ground Floor Living Room Make: Medi-Tec Model: Classic Serial: 2000172843054 S.W.L. 175kg Date of Manufacture 26/10/2016 Installed: Date of Last Examination 14/11/2025 Next LOLER due 14/11/2026</p>
<p>Completion Notes</p> <p>Serviced hoist. Batteries are getting low but still within the 60% threshold. Batteries retained charge after weight test. Checked charge rate. Checked mast, boom and casters. Weight tested to 175kg. Reset service counter. Inspection passed</p> <p>Parts Used</p>	<p style="text-align: center;">Further Work Required</p>

VAT Exemption Certificate - Eligibility declaration by an individual

I declare that i am chronically sick or disabled and suffering from _

I am receiving a qualifying Product or Service/Maintenance visit supplied to me for my domestic or personal use. I claim relief from Value Added Tax under group 12 of schedule 8 of the Value Added Tax Act 1994. The work has been completed.

Customer Name: AIG Control Systems Ltd

Customer Sign:

Date: 14/11/2025

Rep Name: Wayne May

Rep Sign: 

Date: 14/11/2025



LOLER Checks

Customer

AIG Control Systems Ltd
Otter Pool Barn
West Stowford Barton Barn,
West Stowford,
Ilfracombe
Devon, England
EX34 8PJ

C1027

Visit Details

Visit Code: J94541
Order No: TVMC/14/11/2025/Classic
Consultant: AIG Control Systems Ltd
Contact:
Tel:

Make: Medi-Tec
Model: Classic
Serial: 2000172843054

Equipment Examined

Question	Answer
Asset Number if present:	na
What is S.W.L. (Kg)?	175kg

Type of Examination

Question	Answer
Initial Examination	y
Examination after Reinstallation	n
Examination within a 6 -month Interval (9.3.a.i)	y
Examination In accordance with a written scheme	y
Exceptional Examination(Please Give Details)	weight tested to 175kg

Examination Points

Question	Answer
Battery Voltage and Amps (1)- Serial Number of Battery	12.34v 1.9ah
Battery Voltage and Amps (2)- Serial Number of Battery	12.31v 1.9ah
Check Directional Controls/Function	y
Check Charger & Docking Points	26.6v
Emergency Lowering Manual	y
Powered Emergency Lowering	y
Emergency Stop Switch	y
Check Labelling including SWL/ Service Dates	y
Are User Instructions Present?	y
Weight test/Load (kg)	175kg
Pictures taken of weight test?	y
Check Boom to Mast Attachment	y
Check Pump or Actuator Attachment Points	y
Check Mast Security in Base inc. Locking Mechanism	y
Check Alignment of Legs	y
Leg Linkages for Wear/ Security/ Smooth Operation	y
Oil Levels/ Hydraulics	na
Check Carabiner for excessive wear	y
Check Weld Points	y

Castors for Wear/ Security/ Smooth Operation/ Effective Brakes	y
Measurement of castor if required	na
Life cycle of actuator (use calculator to approximate usage)	830
Date of Actuator	26.10.2016
Actuator-full range Wear/ Resistance/Abnormal noise/ Emergency Operation	y
Reset the service light (Molift)	na
P.A.T. Test - Pass/Fail	pass
Dangerous Defects Found?	n
If Yes Please detail items which are safety critical now	n
Potentially dangerous defects (not safety critical) Found?	n
If Yes Please detail for monitoring	n

Consent Request

Question	Answer
Does client give consent for any pictures taken to be used on social media?	na

Rep Name: Wayne May

Rep Sign: 

Date: 14/11/2025







