

## LOLER Report

<p style="text-align: right;"><b>C1027</b></p> <p>AIG Control Systems Ltd Otter Pool Barn West Stowford Barton Barn, West Stowford, Ilfracombe Devon, England EX34 8PJ</p> <p>07970 538998</p>	<p style="text-align: center;"><b>Visit Details</b></p> <p>Visit Type: LOLER Order No: 94540/ Date/Time: 14/11/2025 10:54 Time Completed: 11:48</p>
<p><b>Job Details</b></p> <p>LOLER inspection visit due, 2025</p>	<p style="text-align: center;"><b>Lift Details</b></p> <p>Asset /Location: 1st Floor Landing Make: Medi-Tec Model: Classic Serial: 2000171743054 S.W.L. 175kg Date of Manufacture 26/10/2016 Installed: Date of Last Examination 21/06/2024 Next LOLER due</p>
<p><b>Completion Notes</b></p> <p>Serviced hoist. Batteries are still good and holding charge after weight test. Checked charge rate. Checked mast, boom and casters. Weight tested to 175kg. Reset service light. PAT tested charge lead. Refitted charge lead grommet to control box. Inspection passed</p> <p><b>Parts Used</b></p>	<p style="text-align: center;"><b>Further Work Required</b></p>

### VAT Exemption Certificate - Eligibility declaration by an individual

I declare that i am chronically sick or disabled and suffering from \_

I am receiving a qualifying Product or Service/Maintenance visit supplied to me for my domestic or personal use. I claim relief from Value Added Tax under group 12 of schedule 8 of the Value Added Tax Act 1994. The work has been completed.

**Customer Name:** AIG Control Systems Ltd

**Customer Sign:**

**Date:** 14/11/2025

**Rep Name:**Wayne May

**Date:** 14/11/2025

Rep Sign: *W*



## LOLER Checks

### Customer

**C1027**

AIG Control Systems Ltd  
Otter Pool Barn  
West Stowford Barton Barn,  
West Stowford,  
Ilfracombe  
Devon, England  
EX34 8PJ

### Visit Details

Visit Code: J94540  
Order No: TVMC/14/11/2025/Classic  
Consultant: AIG Control Systems Ltd  
Contact:  
Tel:

Make: Medi-Tec  
Model: Classic  
Serial: 2000171743054

### Equipment Examined

Question	Answer
Asset Number if present:	na
What is S.W.L. (Kg)?	175kg

### Type of Examination

Question	Answer
Initial Examination	y
Examination after Reinstallation	n
Examination within a 6 -month Interval (9.3.a.i)	y
Examination In accordance with a written scheme	y
Exceptional Examination(Please Give Details)	weight tested to 175kg

### Examination Points

Question	Answer
Battery Voltage and Amps (1)- Serial Number of Battery	12.50v 2.6ah
Battery Voltage and Amps (2)- Serial Number of Battery	12.47v 2.7ah
Check Directional Controls/Function	y
Check Charger & Docking Points	26.6v
Emergency Lowering Manual	y
Powered Emergency Lowering	y
Emergency Stop Switch	y
Check Labelling including SWL/ Service Dates	y
Are User Instructions Present?	y
Weight test/Load (kg)	175kg
Pictures taken of weight test?	y
Check Boom to Mast Attachment	y
Check Pump or Actuator Attachment Points	y
Check Mast Security in Base inc. Locking Mechanism	y
Check Alignment of Legs	y
Leg Linkages for Wear/ Security/ Smooth Operation	y
Oil Levels/ Hydraulics	na
Check Carabiner for excessive wear	y
Check Weld Points	y

Castors for Wear/ Security/ Smooth Operation/ Effective Brakes	y
Measurement of castor if required	y
Life cycle of actuator (use calculator to approximate usage)	851
Date of Actuator	22.09.2015
Actuator-full range Wear/ Resistance/Abnormal noise/ Emergency Operation	y
Reset the service light (Molift)	na
P.A.T. Test - Pass/Fail	pass
Dangerous Defects Found?	n
If Yes Please detail items which are safety critical now	n
Potentially dangerous defects (not safety critical) Found?	n
If Yes Please detail for monitoring	n

**Consent Request**

Question	Answer
Does client give consent for any pictures taken to be used on social media?	na

**Rep Name:** Wayne May

**Rep Sign:** 

**Date:** 14/11/2025







